Vivdiam

It's More Than Beautifying, It's An Experience...



- · Cleanses all Dirt, Oils and Makeup Residues
- Reduces Acne Breakouts and Blackheads
- · Reduces Fine Lines and Wrinkles
- Exfoliates Dead Skin Cells for Younger-looking Skin
- Shrinks Large Pores
- Enhances Deep Absorption of Skin Care Products
- Revitalises Dull Skin
- Improves Overall Skin Tone
- Stimulates Collagen Production

PRODUCT SPECIFICATION:

- Power Supply
- Unit Size
- · Unit Weight
- Material of Main Cons
- Material of Main Case
- Vibration Frequency
- Color
- Accessories
- Waterproof Level
- Working Temperature
- Gift Box

- : DC3.7V (200mAx1)
- : L84 x W100 X H30 mm
- : 80a
- : ABS + 100% Food Grade Safe Silicone
- : 7000 rpm
- : Blue, Pink, Purple, Red
- : USB Cord, Manual
- · IPX5
- : Normal Atmospheric Temperature
- : L173 x W115 x H45 mm

PARTNER WITH US!

Enjoy the Benefits of Our "Alliance Program"









Imported, Distributed and Marketed By:

ProExcell International Plt (LLP0011724-LGN)

Pusat Perdagangan Icon City, T1-13A-07, No. 1B, Jalan SS8/39, SS8, 47300 Petaling Jaya, Selangor Darul Ehsan, Malaysia Tel: +60 17-6711 832 / +60 18-3718 311

Web: vivdia.proexcell.com.my



Vivdia

Multi-Function
Beauty Facial Massager



YOUR ULTIMATE PROBEAUTIFYING SYSTEM

Get Youthful, Healthier, Radiant Skin

EXPERIENCE THE CHANGE TODAY

600 Smooth Massage Contacts In the Top
 Use it to tackle those difficult to clean areas, especially
 around the nose, ears, and around the hairline with a
 relatively strong vibration and a little tingling sensation
 during usage.

 More than 1400 Soft and Smooth Clean Contacts for Skincare Applications

Feel the softness of the contacts on your face, the contact vibration massage is suitable for sensitive skin. Uses thousands of pulsations per minute to penetrate your skin effectively.

• 15 Different Speed Vibration Modes

15 different levels of personalized adjustable massage intensity for all skin types.

skin typ

• Eco-Friendly
USB rechargeable,
quick charge for 30-60
minutes. Turns off
automatically after 15
minutes of inactivity.

 No Replacement Brushes Needed Never replace another brush head again,

Never replace another brush head again, and brushes continue to stay hygienic.

Certified Waterproof (IPX5)
 Standard

The whole fuselage can be washed, and includes a small safety cover over the charge port, is waterproof (IPX5), safe and convenient to use.

Provides Comprehensive Cleansing

Dispel impurities and oil in difficult-to-clean corners with ease.

Aids Skin Cells Absorption Within Seconds

Aids absorption into skin cells and stimulates nutrient and collagen production within seconds of all skin care ingredients after application.

100% FDA-certified Food Grade Safe

In case of ingestion, the product will not cause any harm to the person.

Global Assurance

The product is EMC, FCC, CE and ROHS certified.

How to Use?

1. Cleaning and Cleansing Mode











Step 1: Smear cleanser

Step 2: anser Clean both cheeks

ep 2: Ste ean both Cle eeks

Step 3: Step 4: Clean forehead Clean nose

Clean following parts of the

2. Firming, Whitening and Soothing Mode









3

Between the eyebrows

Right eye tail

Step 3: Right nose lips


Step 4: Left nose lips


Step 5: Left eye tail

Usage Caution

- Please don't use this product for the sensitive area around the eyes.
- Please don't use this product more than twice a day.
- Please do not use hot water (up to 40 degrees Celsius to clean the product).
- If there is broken skin, gaping wound, skin disease or severe acne, sunburn, skin infections, do not use this product.
- If you have sensitive skin, please pay attention to the massage intensity.
- We recommend that you don't clean your skin for over 1 minute at a time, in case of excessive cleansing.
- Directly rinse with clean water after cleansing is completed.

How to Charge?





- 2 indicator light flashing during charging.
- The indicator light will be constantly on once it has been fully charged.
- · Quick charging time: 30 60 minutes.
- · 3-4 hours usage time after fully charged.